
CELLINA S
AGELESS IMAGINATION
Dual fractional RF system
Precise dual wave technology

CELLINA S

The evolutionary entity in RF combination


Since its initial launch in 2014, our efforts in technological advancements for our RF devices have evolved dramatically over the past decade. We proudly present Cellina S, the 3rd generation RF device that offers more indications, comfortable procedures, and faster recovery in daily life.

Advanced RF-microneedling combined with fractional RF provides optimal customized effects for skin rejuvenation, lifting & tightening, skin tone improvement, and acne treatment.




Ablative Tx






Gentle Pro treatment penetrates the skin with tiny microneedles, precisely delivering controlled RF energy to the deeper layers of the skin.


Non-Ablative Tx





Subtle treatment stimulates new collagen growth and resurfaces the skin by delivering bipolar RF energy to the skin via multi-electrode pins.


Disposables



∗ Option

SAP Technology


Safe Anti-Pain Technology minimizes pain by delivering customized, stable and ideal RF energy deep into the skin and achieves optimal results for each treatment.



TECHNOLOGY

User-friendly custom solution with advanced design

Insulated and non-insulated microneedles come in two types, allowing for tailored procedures without skin damage based on skin type and individual patients. Specially designed, incredibly thin gold-coated needles, along with upgraded motor, facilitate painless and fast procedures, as well as quick recovery in daily activities.





Precise and selective application

The dual wave combination of Cellina S is capable of implementing various treatment modes to the maximum extent.

Continuous Wave : Irradiating the dermis layer with continuous energy generates high heat instantly, causing dermal layer coagulation and maximizing skin lifting and tightening effects.

Pulse Wave : Various pulse modes deliver energy to the target area without damaging the surrounding tissue with minimal pain, effective for treatments such as pigmentation, acne, and melasma.





Ageless Imagination

The treatment was performed without causing any discomfort or intense pain, making it an enjoyable experience.

- Elsa Lindberg, Sweden

The treatment provided me with a noticeable improvement in the texture and tone of my skin.

- Haruka Tanaka, Japan

I used to worry a lot about my neck wrinkles, but now I don't care anymore after getting the treatment.

- Sophie Dubois, France



PRODUCT SPECIFICATIONS

Product Characteristics	Specifications
Technology	Radiofrequency
Frequency	1 Mhz / 2 Mhz
Power	max. 60W
Total RF treatment levels	up to 10 levels
Pulse intervals	Continuous / Pulse
Repetition rate	0.5~3 sec (0.5 sec step)
Display	Touch LCD Display
System cooling	air-cooled
Weight	30 kg
Dimensions (including arm)	390mm(W) x 520mm(L) x 1400mm(H)

DAOL MED
Office: +82 (0)2 6925 2199, Fax: +82 (0)2 6925 2450
Mail: sales@daolmed.co.kr www.daolmed.co.kr

WAYWELL
WAY TO WELLNESS
Office: +82 (0)2 6925 0122, Fax: +82 (0)505 065 0122
Mail: info@waywellgroup.com www.waywellgroup.com

CE 0068